


SCORPION SUTURE PASSER

- Scorpion suture passer is used in shoulder arthroscopy for rotator cuff and labral repair.
- Designed for simplicity, the easy-to-use Scorpion side-loading suture passer supports more successful suture passes.
- The front side of the bottom jaw is solid to prevent suture and soft tissue from entering the bottom jaw when positioning onto the tendon.
- The suture-retrieving top jaw helps seamlessly pass and retrieve #2 Fiber suture.
- The Scorpion suture passer has a quick-loading bottom jaw and suture-retrieving top jaw.
- This intuitive suture passer's streamlined, low-profile jaws are designed for easily maneuvering the glenohumeral space, fit through a 7.0 mm cannula, and are ideal for passing and retrieving #2 Fiber suture through a single portal.



Order Informations

Ref Code SCSP1000 SCORPION SUTURE PASSER LONG
SCMP1001 SCORPION SUTURE PASSER SHORT

NeoSys

www.neosysmed.com



SHOULDER
TENDON AND
LIGAMENT REINSERTION