


FIBERTAPE SUTURE

-> FiberTape suture offers broad compression and increased resistance to tissue pull-through, making it the leading choice for knotless repairs.
-> This ultra-high-strength suture has a long-chain polyethylene structure similar to Fiber suture.
-> Its broad footprint is ideal for repairs to degenerative cuff tissue where tissue pull-through may be a concern.
-> FiberTape suture is 30% more resistant to tissue pull-through compared to round #2 suture.



Order Informations

Ref Code FBT1490 - FIBER TAPE 1.4MMX90CM
FBT2090 - FIBER TAPE 2.0MMX90CM
FBT4090 - FIBER TAPE 4.0MMX90CM

NeoSys

www.neosysmed.com



SHOULDER
TENDON AND
LIGAMENT REINSERTION