


POLYMER CERCLAGE CABLE SYSTEM

- The cable is made from a high-resistance polymer.
- The cable can be used with a special instrument set, taking shape as desired.
- There are 2 scales and 1 torquing device in the instrument set.
- Finds clips made of titanium that allow to lock on the cable.
- The cable is color-coded and has a characteristic that the light can absorb and shine.
- The cable is ready for use in the sterile box.
- The cable length: 110 cm, dia: 1.6 mm.
- It can be passed through all cable plates or can be used alone.



Order Informations

Ref Code KBL.ORTH.001 POLYMER CERCLAGE CABLE

NeoSys

www.neosysmed.com

