


# SELF-PUNCHING KNOTLESS ANCHOR

- The self-punching 4.75 mm knotless anchor combines a titanium tip with a PEEK and titanium anchor body to eliminate the need for prepunching a bone socket.
- This self-punching design can help save valuable OR time, while increasing the precision of the final construct.
- The self-punching anchor can be combined with FiberTape suture to complete a SpeedFix repair or SpeedBridge knotless rotator cuff repair.



## Order Informations

Ref Code 16645185 KNOTLESS ANCHOR PEEK+TI 4.75MM

# NeoSys

[www.neosysmed.com](http://www.neosysmed.com)



**SHOULDER**  
TENDON AND  
LIGAMENT REINSERTION