

# FEMORAL FIXATION BUTTON CL

→ Pre-measured loop size button is made of:

- ... 1 titanium button (dimensions 12x4 mm)
- ... Polyester PET braided pre-measured loop

→ 2 pre mounted UHMWPE sutures (green and white)

## → Indications

- ... Tendon or ligament fixation in reconstruction surgeries especially for ACL repair.

## → Loop sizes

- ... 6 sizes are available and meet requirements of conventional ligamentoplasty surgeries (**15, 20, 25, 30, 35, 40mm**)

→ This button implant, which is clinically proven in effectiveness and reliability, is designed by orthopedic surgeons for anterior cruciate ligament surgery, femoral tunnel hamstring graft cortex hanging detection.

→ Pulling sutures are different colors to avoid confusion during processing.



KNEE

## Order Informations

Ref Code FKL.ART15 - FEMORAL FIXATION BUTTON CL 15MM


Ref Code FKL.ART20- FEMORAL FIXATION BUTTON CL 20MM

Ref Code FKL.ART25 - FEMORAL FIXATION BUTTON CL 25MM

Ref Code FKL.ART30 -FEMORAL FIXATION BUTTON CL 30MM

Ref Code FKL.ART35 - FEMORAL FIXATION BUTTON CL 35MM

Ref Code FKL.ART40 - FEMORAL FIXATION BUTTON CL 40MM



# NeoSys

[www.neosysmed.com](http://www.neosysmed.com)