

PEEK SUTURE ANCHORS

- NEOSYS PEEK anchor is a fully threaded suture anchor featuring dual threads to maximize cortical and cancellous fixation, has a rounded tip to protect the sutures and to facilitate the insertion.
- Particularly suitable for repairing rotator cuff and associated pathologies, this anchor is available with two or three sutures and/or tapes.
- Double cancellous and cortical thread.
- Maximize fixation strength by increasing the interaction between cortical and cancellous bone.
- Greater torsional torque when inserted.
- Maximum tensile strength, therefore pull-out.
- Optimal internal footprint of the screwdriver.
- More important contact surface between the anchor and the screwdriver for better torque and guided insertion.
- Sizes:**
3.5 mm, 5.0 mm
- Suture Configurations:**
1 suture or 2 sutures
- Adaptable:**
Suitable for multiple surgical techniques.

It is specially designed for the fixation of soft or prosthetic tissues to the bone.

- NEOSYS PEEK Suture Anchors provide surgeons with pullout strength, high-strength Fiber Suture, and multiple hole preparation options for any bone quality.
- Transitional thread design enables secure fixation with distal “cutting” and proximal “locking” threads, which means easy insertion and stronger pullout strength even in poor quality bone.
- Loaded with fiber sutures inside the body of the anchor means stronger knots.
- General hole preparation instrument offers solutions suitable for all qualities of bone.
- With a modulus similar to cortical bone, VestaPeek® polymer from Evonik® exhibits an ideal combination of strength, stiffness, and toughness, together with biocompatibility.



Order Informations

Ref Code 16450017 PEEK SUTURE ANCHOR 5.0MM WITH NEEDLE

Ref Code 16435017 PEEK SUTURE ANCHOR 3.5MM WITH NEEDLE

NeoSys

www.neosysmed.com



SHOULDER
TENDON AND
LIGAMENT REINSERTION