


DOUBLE BUTTON FOR SYNDESMOSIS REPAIR

- **SYNDESMOSIS Double Button** with Adjustable Loop is designed for syndesmosis.
- It's an important innovation that provides great convenience for surgery in syndesmosis operations, which have a very important place in orthopedic surgical indications and that satisfies a desired fixation at the maximum level.
- Additionally, it's used for syndesmosis according to surgical indications in important joints such as the ankle, knee, and wrist.
- Double Button is a precontoured, titanium button that allows the use of multiple Fiber Sutures for syndesmotic reduction, providing a construct that is twice as strong as existing repair devices.
- Since the buttons are attached to the Fiber sutures independently, only suture material is passed through the fibula and tibia tunnels, allowing the repair to be secured with minimal disruption to the surrounding bone structure.



Order Informations

Ref Code 8626T

NeoSys

www.neosysmed.com