


GRAFT PASSING GUIDE PIN

- The 2.4mm Graft Passing Guide Pin is specifically engineered for fast, accurate, and minimally invasive referencing of the Aperture to Cortex (AC) length.
- The fluted drilling head and tapered shaft allow for a tactile feel in determining the AC length.
- Surgeons can confidently confirm the femoral tunnel length in the early stages of the reconstruction.
- This pin helps save time and minimizes potential complications.



Order Informations

Ref Code 1297L - 2.4MM PASSING PIN

NeoSys

www.neosysmed.com