

ALL SUTURE SOFT ANCHORS

It is specially designed for the fixation of soft or prosthetic tissues to the bone.

→ SOFT ANCHOR are impacted anchor of small diameter exclusively composed of textile fibers.

→ The polyester body of the anchor is associated with one or two UHMWPE suture threads which are mounted on a disposable screw-driver allowing insertion into the bone.

→ The traction on the suture threads makes it possible to modify the shape of the body of the anchor which thus locks itself behind the cortical bone and ensures an effective anchorage.

→ **All Suture soft Anchor - 1.8 , 2.9 mm:**

Most commonly utilized for soft tissue to bone fixation in the shoulder. Designed with a custom braid polyester sleeve and double loaded with #2 UHMWPE - ULTRAFIRE FIBER suture.

→ **All Suture soft Anchor - 1.0 - 1.5 mm:**

Most commonly utilized for soft tissue to bone fixation in the shoulder and small joints. Designed with a custom braid polyester sleeve and single loaded with #2 UHMWPE - ULTRAFIRE FIBER suture.



Order Informations

Ref Code 9901002 Soft Anchor 1.5mm

Ref Code 9901004 Soft Anchor 1.8mm

Ref Code 9901005 Soft Anchor 2.9mm

Ref Code 9901001 Soft Anchor 1.0mm



HAND
INSTABILITY



SHOULDER
TENDON AND
LIGAMENT REINSERTION



ANKLE
REINSERTION

NeoSys

www.neosysmed.com