


ADJUSTABLE LOOP BUTTON

- The adjustable femoral cortical fixation system provides adjustable cortical fixation for anterior cruciate ligament reconstruction.
- It offers strong pulling force.
- Lift button removes the need for multiple implant sizes.
- The lift button Femoral Cortical Fixation System is designed for graft femoral cortical attachment.
- The implant has a variable knot length structure that allows the tendon within the tunnel to be adjusted with optimum tension.
- Single implant option for entire tunnel length UHMWPE sutures for pulling.



Order Informations

Ref Code AR1588T

NeoSys

www.neosysmed.com