

SWIVELOCK ANCHORS

- Fully threaded knotless SwiveLock anchors are designed for use with SutureTape, Fiber Sutures, Fiber Tapes, and soft-tissue grafts in repair and reconstruction techniques.
- Tension is visualized, adjusted, and locked into position with the SwiveLock anchor body.
- Available in PEEK material, SwiveLock anchors are cannulated and vented to minimize material.
- 2.9mm SwiveLock anchor is available for instability repair and provides very high pullout and insertion strength while saving time.



Order Informations

- Ref Code SWL290 2.9mm swivelock anchor
- Ref Code SWL0350 3.5mm swivelock anchor
- Ref Code SWL0450 4.5mm swivelock anchor
- Ref Code SWL0550 5.5mm swivelock anchor

NeoSys

www.neosysmed.com



KNEE



HIP



HAND



SHOULDER



ANKLE