

MINI SUTURE ANCHOR (2.0 MM)

- 2.0mm Mini suture anchor is designed for efficiency and ease of use, thereby optimizing patient outcomes.
- Our Mini Titanium Screw-In Suture Anchors are recommended to address a variety of small-joint applications due to their small diameter and tremendous pull-out strength.

Key Features & Benefits:

- Threaded design provides optimal cortical fixation and reduces the risk of “pull-back.”
- Made of 6Al-4V ELI titanium alloy, providing the strength and stiffness required for successful application.
- Self-tapping threads.



Order Informations

Ref Code 16120015

NeoSys

www.neosysmed.com



HAND



ANKLE