

AC DOUBLE BUTTON SYSTEM

- The AC double button system introduces a novel technique for the reconstruction of complete acromioclavicular joint separation.
- It utilizes buttons with CL material, which has been modified for use in the shoulder.
- The device is placed through holes in the coracoid and clavicle, reproducing the course of the conoid portion of the coracoclavicular ligament.
- The buttons with CL material provide strength and stiffness exceeding that of the native anatomy, ensuring a stable reduction.
- The procedure is simple, has low morbidity, and can be easily adapted to an arthroscopic technique.



Order Informations

Ref Code 2257

NeoSys

www.neosysmed.com



SHOULDER
TENDON AND
LIGAMENT REINSERTION